

## AQUATRAM™ Pool / Spa Verification Guide

Thank you for considering the Pentair AquaTram<sup>™</sup>, PN 11200, ADA compliant pool lift. Please take the time to verify the below information to ensure that the AquaTRAM is appropriate for your facility *prior to ordering*. NOTE: It is the responsibility of the buyer to confirm both that this lift is appropriate for your pool deck and gutter profile; and that the resulting installation is ADA compliant. For assistance call 1-800-831-7133.



## Step 1.


Is your pool profile similar to one of the four shown?

☐ YES ☐ NO\*


Is your deck-to-water 9" or less?

☐ YES ☐ NO\*

<sup>\*</sup> If your deck-to-water is greater than 9" but you answered "YES" to the other questions, ask about our AquaTRAM **Deep Draft** model.



Caution: The AquaTram lift is not designed for pool profiles with curbs, walls or parapets. However, it <u>may</u> be able to clear very small curbs (up to 6" wide x 4" high) depending on curb location relative to the pool wall. Please do not order the AquaTRAM for use on a pool with even a small curb without confirmation of clearance from Pentair Applications Assistance: 1-800-831-7133



<u>IMPORTANT!</u> For assistance in understanding the ADA please visit <u>www.access-board.gov/recreation/guides/pools.htm</u> or contact Pentair by calling 1-800-831-7133. Should a lift need to be returned due to improper selection: All returns are subject to our written approval and must be accompanied by "Return Goods Authorization" form or number issued by us prior to returning the lift. Unauthorized returns will not be accepted. Freight must be prepaid. Do not request authorization to return material or credit which has been in your inventory over sixty (60) days and which is not in original resalable condition. Lifts authorized to be returned for credit are subject to a minimum 15% handling charge.