

# HAFFMANS CPT

## CO<sub>2</sub> PURITY TESTER

### GENERAL PRODUCT INFORMATION


In the beer and beverage industries, the quality of carbon dioxide (CO<sub>2</sub>) is a determining factor in the product's quality and taste. A low oxygen (O<sub>2</sub>) content is crucial to avoid oxygen pickup in beverages. The CO<sub>2</sub> collected from fermentation is liquified in order to reduce the space needed for storage and to purify the CO<sub>2</sub> from non-condensable gases - O<sub>2</sub> and nitrogen (N<sub>2</sub>).

The CO<sub>2</sub> Purity Tester, type CPT, is used for routine controls of the air content of commercial CO<sub>2</sub> gas and CO<sub>2</sub> gas from fermentation. The principle is based on CO<sub>2</sub> absorption by caustic.

The measuring burette of the CPT is properly flushed with CO<sub>2</sub> gas. The CO<sub>2</sub> sample is separated in the measuring burette and the connection is made with the caustic burette. The CO<sub>2</sub> is completely absorbed by the caustic and the remaining gas volume is read on the scale of the measuring burette after turning the instrument on its backside.

Two executions are available:

- Type CPT 99-100, measuring burette scale from 99 - 100% v/v, for measuring the end purity of commercial CO<sub>2</sub> gas or CO<sub>2</sub> gas from fermentation after purification
- Type CPT 50-100, measuring burette scale from 50 - 100% v/v, for measuring the intake purity of CO<sub>2</sub> gas from fermentation.



### CUSTOMER BENEFITS

- **Perfect product control.**
- **Cost saving:**
  - less labor intensive - multiple measurements without the need to remove caustic.

### APPLICATIONS

- **Quality control:**
  - of gas distribution systems.
  - of commercial gas.
  - of incoming CO<sub>2</sub> gas from fermentation.
  - of CO<sub>2</sub> recovered and purified from fermentation.
  - in front of the carbonator to monitor the gas quality dosed.
- **Optimizing CO<sub>2</sub> recovery plants:**
  - behind the activated carbon filter/dryer.
  - in the purge from the CO<sub>2</sub> liquefier to control the blow-off quantity automatically, preventing CO<sub>2</sub> loss as well as high concentrations of O<sub>2</sub> from dissolving CO<sub>2</sub>.
  - behind the CO<sub>2</sub> evaporator to monitor O<sub>2</sub> in the CO<sub>2</sub> distribution system.

# HAFFMANS CPT

## CO<sub>2</sub> PURITY TESTER

### TECHNICAL DATA

#### CPT 99 - 100%

##### Measuring range

99.0 - 100.0% v/v CO<sub>2</sub>

##### Graduated scale

99 - 100% 0.02% v/v foreign gas

##### Accuracy

99 - 100% 0.01%

#### CPT 50 - 100%

##### Measuring range

50.0 - 100.0% v/v CO<sub>2</sub>

##### Graduated scale

94 - 99.8% 0.1% v/v foreign gas

55 - 85% 2.0% v/v foreign gas

##### Accuracy

94 - 99.8% 0.05%

55 - 85% 1.0%

##### Volume

170 ml caustic solution

##### Caustic concentration

30% (m/m) or (wgt/wgt) or g/100g KOH

with approx. 2% ethanol

20% (m/m) or (wgt/wgt) or g/100g NaOH

with approx. 2% ethanol

##### Dimensions

430 x 110 x 330 (L x W x H mm)

##### Weight

Ca. 1 kg

### SCOPE OF SUPPLY

- CO<sub>2</sub> Purity Tester\*
  - Type CPT 99-100 or
  - Type CPT 50-100
- Instruction manual

\*When ordering, please state which instrument you require.

### OPTIONS

- Certificate of measurement



#### HAFFMANS B.V.

P.O. BOX 3150, 5902 RD VENLO, NETHERLANDS FOODANDBEVERAGE.PENTAIR.COM

All Pentair trademarks and logos are owned by Pentair. All other brand or product names are trademarks or registered marks of their respective owners.

Because we are continuously improving our products and services, Pentair reserves the right to change specifications without prior notice.

Pentair is an equal opportunity employer.

©2026 Pentair, All Rights Reserved.

